

Question:

Do you have some instructions for how to attach the Velcro Limb Straps? I think we need a diagram! I don't know if velcro side touches the skin, but it seems like the smooth side should, but I can't figure out how to get the thing to stay put....

Does the metal disk need to touch the skin??

BioTrode™ Limb Straps and Wires

The rubber side goes against the skin.

Start with the end that does not have the magnet. Wrap around and you may be able to wrap twice, overlapping, which will give you more skin contact. Then the small Velcro® tab at the other end simply attaches to the soft Velcro® and you are strapped up.

Next, attach the magnets on the wire (color coded) to the magnet on the strap.

No metal is ever to touch the skin in order to prevent EMF burns.



Go to the website and take a look www.BioTrodeTech.com



Be sure to take a look at the **FAQs** tab while you are there and:

On the homepage you can sign up for news letters, promotional coupons (We don't send much out so you won't be flooded).

Sorry for the confusion. We are currently making a DVD on how to accessorize but is slow going.

If you have further questions, call me.

Joni

(425) 277-5121

Reply:

These instructions are PERFECT! I didn't see this on your website, but I might not have looked in the right place. I think if this was typed on a slip in with the packing papers, it would be wonderful :) Oh these are so much more comfy than the regular ones that come with the device! They feel much more secure and I love that they unsnap for the dreaded bathroom break....